HIRSCHMANN®

Product data sheet

Sealing Wax Plates

Manufacturer	HIRSCHMANN LABORGERÄTE HAUPTSTRAßE 7-15 • 74246 EBERSTADT
Item number	9120102
Item description	SEALING WAX PLATES
Product image	
Product features	• Numbered from 1 to 24 • 24 storage options • For sealing capillaries, e.g. hematocrit capillaries
Packaging units	1 pack of 10 pieces
Packaging	1 pack is contained in a white product box 
Area of application	The plates are used as a sealing method in all applications where capillaries are used, for example in hematocrit capillaries.
Compatibility	The sealing wax plates can be used for all Hirschmann® capillaries.
Dimensions (mm) and weight (g) of a plate	Length: approx. 118 mm Width: approx. 42 mm Height: approx. 7 mm Weight: approx. 18.4 g
Materials	Plate: PP Filling material: plasticine
Sterility	The product is not sterile.
Usage period & Storage	We recommend a usage period of approximately 3 years from the date of purchase. If stored correctly between 10–25 °C, the product can be used for longer.

DOMINIQUE DUTSCHER SAS

Germany
Hirschmann Laborgeräte GmbH & Co. KG
Hauptstraße 7-15
74246 Eberstadt
Phone +49 7134 511-0

For more information, please visit:
www.hirschmannlab.de

USA
Hirschmann, Inc.
Laboratory Equipment
2700 Holloway Road, Suite 104
Bluegrass Industrial Park
KY 40299 Louisville
Phone: +1 502 2400287

DOMINIQUE DUTSCHER SAS

The information presented in the above data sheet represents standard values. Hirschmann assumes no responsibility for their accuracy. Users are solely responsible for the suitability of the products, the conduct of preliminary investigations, and compliance with legal regulations. Use of the products for implantations or injections is not permitted. Technical specifications may be changed without prior notice; Hirschmann reserves this right.
©Hirschmann Laborgeräte GmbH & Co. KG