

Staubschutzkappe GL45
Zeichnung: 0817201-GZA-0GL45-..

Gewinde:
DIN 168-GL 45x4; 3 Gang
Zeichnungs: 0303803-GZA-01694-..





100 ml ±5%
DURAN
Made in Germany
Retrace Code
00000000

DURAN PURE 80
60
40

Farbe: weiß
0388101-SZ1-A

(22)

DOMINIQUE DUTSCHER SAS

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.
The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

A	DWK-Logo		21-AUG-2017	Weindorf
Index	Aenderung / changes		Datum / date	Name

zul. Abweichung allowed dev.		
Oberfläche surface		
	Datum / date	Name
bearb. drawn	13.11.2012	Weindorf
geprüft checked	31-08.2017	Fries

DURAN WHEATON KIMBLE
DURAN Produktions GmbH & Co. KG
www.DWK-LifeSciences.com

Werkstoff / material	Stück quantity	Gewicht weight	Pos. position
DRAN			

Benennung / title
Laborflasche 100 ml
DURAN PURE
Art.-Nr. 81 801 24 06

Maßstab scale

Form-Nr. mould number: 1694/100 — A3

Zeichnungsnummer / drawing number: 0340219-GZA-01694 — Revision: A — Blatt/sheet: 1 /1

Ers. / repl.f: | Ers. / repl.d:

Schutzvermerk ISO 16016 beachten. / Refer to protection notice ISO 16016. — INV