

FOR CUSTOMER REFERENCE ONLY


THIS PRODUCT CONFORMS TO:
ANSI/SBS 1-2004 FOOTPRINT DIMENSIONS
ANSI/SBS 4-2004 WELL POSITIONS



PRODUCTS

PART #	PLATES / KIT	DESCRIPTION
4483354	20	MicroAmp® EnduraPlate™ Optical 96-Well Clear
4483343	20	MicroAmp® EnduraPlate™ Optical 96-Well Blue
4483349	20	MicroAmp® EnduraPlate™ Optical 96-Well Green
4483350	20	MicroAmp® EnduraPlate™ Optical 96-Well Red
4483395	20	MicroAmp® EnduraPlate™ Optical 96-Well Yellow
4483355	5	MicroAmp® EnduraPlate™ Optical 96-Well Multicolor
4483352	500	MicroAmp® EnduraPlate™ Optical 96-Well Clear
4483356	500	MicroAmp® EnduraPlate™ Optical 96-Well Multicolor





UNLESS OTHERWISE SPECIFIED: DIMENSIONS AND TOLERANCES ARE IN MILLIMETERS AND PER ASME Y14.5.
TOLERANCES ARE:
ANGLES ±.5°
X. ±.5
X.X ±.25
X.XX ±.13
X.XXX ±.05
0.8
REFERENCED 3D FILE REVISION: A-09
AGILE NUMBER:
METRIC

MATERIAL: WELLS: POLYPROPYLENE SHELL: POLYCARBONATE
FINISH: SEE NOTE 8
DRAWN: S. JAFFE DATE: 30 NOV 11
CHECKED: CS LIM DATE: 21 MAY 12
QUALITY: DATE:
DOC. CONTROL: DATE:

DESCRIPTION: MICROPLATE, 96 WELL, ENDURAPLATE
SIZE: B SCALE: 1:1 PART NUMBER: 4473712 REV: A-09
DO NOT SCALE DRAWING 3RD ANGLE SHEET 3 OF 3

COMMENTS:
© 2011 Life Technologies Corporation. All Rights Reserved. These drawings and specifications are the property of Life Technologies Corporation and are not to be reproduced, copied, or used as the basis for manufacture or sale of, or any unauthorized purposes related to, the matter depicted herein, without the prior written consent of Life Technologies Corporation.

DOMINIQUE DUTSCHER