Tissue Culture Flask 25 – 300 cm² with Peel-Off Foil



DOMINIQUE DUTSCHER SAS

Features of the **Tissue Culture Flask with Peel-Off Foil:**

- Available in all sizes: 25 - 300 cm²
- Large pinched flap for a firm grip
- Foil is peeled off before access to the cells
- Re-closure of the flask is not possible
- Stacking ring enables stable, non-slip stacking of several flasks
- Optionally available with barrier

Flask is available with filter screw cap only. If a VENT screw cap is used, a possible pressure build-up of more than 0.03 bar inside the flask could press to open the foil.

Type	Product-No.	Growth Area cm²	Version	Lid Opening l x w mm	Volume max. ml	Volume recom. ml	Dimensions l x w x h mm	Material	Qty/ Bag pcs	Qty/ Case pcs
	90028	25	peel-off	38 x 45	15	3 - 8	92 x 51 x 31	PS	5	70
	90078	75	peel-off	87 x 92	65	8 - 22	155 x 87 x 45	PS	2	32
	90153	150	peel-off	100 x 128	100	15 - 45	210 x 122 x 50	PS	3	18
	90303	300	peel-off	140 x 185	200	30 - 85	275 x 170 x 55	PS	2	8
	90653	115	peel-off, barrier	100 x 128	100	15 - 33	210 x 122 x 50	PS	3	18





1. Large flap for easy peel-off
2. Barrier optionally available