Tissue Culture Flask with Re-Closable Lid 115 - 150 cm²






Free access to the cell chamber

DOMINIQUE DUTSCHER SAS

Features of the **Tissue Culture Flask with Re-Closable Lid**:

- A multi-use re-closable liquid-tight lid
- Complete access to flask chamber from above
- Safety edge enables stable, non-slip stacking of several flasks also in combination with the standard TPP tissue culture flasks
- Optionally available with barrier which reduces the growth area to 115 cm²

Flask is available with filter screw cap only. If a VENT screw cap is used, a possible pressure buildup of more than 0.03 bar inside the flask could press open the lid.



Opening and closing mechanism

1 Press here
2 Lift lid

Type	Product-No.	Growth Area cm²	Version	Lid Opening l x w mm	Volume max. ml	Volume recom. ml	Dimensions l x w x h mm	Material	Qty/ Bag pcs	Qty/ Case pcs
	90552	150	lid	105 x 105	100	15 - 45	210 x 122 x 50	PS	3	18
	90652	115	lid, barrier	105 x 105	100	15 - 33	210 x 122 x 50	PS	3	18





1. Secure stacking of multiple flasks
2. Opening and re-closing mechanism