Vacuum Filtration "rapid"-Filtermax 150 – 1000 ml



Quadratic vacuum filtration, conveniently excellent!





1. Filter top and flask are available separately
2. Flask can be used as medium reservoir

Features of the quadratic **vacuum filtration "rapid"-Filtermax:**

- Large quadratic filter surface of 49 cm² and 69 cm² respectively
- Special lamella grid support and premium PES filter membrane, pore size of 0.22 µm with low protein binding capacity resulting in an extreme high flow rate
- Square thus secure footprint eliminates the risk of encompassment of the filter unit by the vacuum tubing
- Space saving medium reservoir
- Tapered form at the neck enables large and small hands a firm and slip-free grip of the "rapid"-Filtermax system
- Each filter set comes with an adapter for the vacuum connection hose with an inside diameter of Ø 6 - 10 mm and a separately packed screw cap

Full flexibility: Choice of filter set, bottle top filter only or flask only, all provided with a GL-45 compatible screw thread.

TPP recommends filling the liquid into the filter first and then applying the vacuum.

Type	Product-No.	Filtration Area	Volume	Pore Size	Dimensions l x w x h	Material	Qty/ Bag	Qty/ Case
		cm²	ml	µm	mm		pcs	pcs
	99150	49	150	0.22	93 x 93 x 103	PS / PES	1	18
	99155	49	150	0.22	90 x 89 x 57	PS / PES	1	36
	99157	-	150	-	93 x 93 x 63	PS	1	24
	99250	49	250	0.22	93 x 93 x 145	PS / PES	1	12
	99255	49	250	0.22	90 x 89 x 75	PS / PES	1	24
	99257	-	250	-	93 x 93 x 87	PS	1	24
	99500	49	500	0.22	93 x 93 x 214	PS / PES	1	10
	99505	49	500	0.22	90 x 89 x 110	PS / PES	1	21
	99507	-	500	-	93 x 93 x 121	PS	1	36
	99560	49	500	0.1	93 x 93 x 214	PS / PES	1	10
	99950	69	1000	0.22	111 x 111 x 285	PS / PES	1	9
	99955	69	1000	0.22	108 x 108 x 143	PS / PES	1	12
	99957	-	1000	-	111 x 111 x 160	PS	1	20
	99960	69	1000	0.1	111 x 111 x 285	PS / PES	1	9

DOMINIQUE DUTSCHER SAS