
Cell Scraper and Cell Spatula



100 % cell harvest

Features of the **Cell Scraper**:

- Scraper heads made from a flexible polyethylene material
- Pivoting scraper head for easy maneuvering
- Three different blade widths
- Small raised knobs on handle for secure grip
- Single-wrapped in blister packaging
- User-friendly dispenser case

Features of the **Cell Spatula**:

- Blade with sharp edges
- Small raised knobs on handle for secure grip
- Single-wrapped in blister packaging
- User-friendly dispenser case

The pivoting scraper head of the cell scrapper reaches every corner of the cell vessel. Pivot the scraper head by slight pressure on the handle with the forefinger thus pushing the handle downwards. A slight move of the handle will twist the head to the required position.

The cell spatula features a specially formed blade with a sharp edge. It is commonly used with the TPP tissue culture flask with re-closable lid or peel-off foil or with large tissue culture dishes.

DOMINIQUE DUTSCHER SAS



Rotation to the required position



Different blade widths available




Cell scraper in a large TPP tissue culture dish

Type	Product-No.	Length	Version	Blade Width	Material	Qty/ Bag	Qty/ Case
		mm		mm		pcs	pcs
	99002	240	scraper "S", rotating	13	PE	1	150
	99003	300	scraper "M", rotating	20	PE	1	150
	99004	362	scraper "L", rotating	25	PE	1	100
	99010	195	spatula, fix	14	PP	1	150

DOMINIQUE DUTSCHER SAS





1. Knobs for secure grip
2. Soft scraper head to protect the cells