PanReac AppliChem ITW Reagents

PRODUCT CODE: 256991

Harris Hematoxylin modified solution (CE-IVD) for clinical diagnostics

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples., In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 00 00 00

PHYSICAL DATA: liquid, garnet, Miscible with water

BIBLIOGRAPHY: Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore (1981) ; H. J. Conn , "Biological Stains", 9th Ed., Williamsand Wilkins , Baltimore , 472 (1977) ; R. García del Moral,"Laboratorio de Anatomía Patológica", 1ª Ed., McGraw-Hill , 129 (1993); E.B. Prosphet , B. Mills , J. B. Arrington , L. H. Sobin ,"Laboratory Methods in Histotechnology", Armed Forces Institute ofPathology , 56 (1992) ; S. Frankel , S. Reitman , A. C. Sonnenwirth ,"Gradwohl#s clinical l •

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,095 kg 1kg~0,913 l

OBSERVATIONS: Filter any sediment present in the solution. •

SPECIFICATIONS:

Suitability for hematoxylin - eosin stain	passes test

DOMINIQUE DUTSCHER SAS