OWIBIO
By Ozyme & Clearline



ELECTROPHORESIS
Owi ExactLadder® DNA 100bp Ready to Load


Reliability. Accuracy. Sucess.

We are with you to improve your research.

DOMINIQUE DUTSCHER SAS

PCR
Polymerase (Standard, Hot Start, High Fidelity), dNTP mix

RT-qPCR
Master Mix, dNTP mix




WESTERN BLOT
ECL, Markers

ELECTROPHORESIS

- **Owi ExactLadder® DNA 100bp Ready to Load**
- Owi Agarose Molecular Grade - Low EEO
- Owi Agarose Standard
- Owi EmeralDye Safe DNA Stain

OWIBIO
By Ozyme & Clearline

Owi ExactLadder® DNA 100bp Ready to Load



5 µl loading (500 ng)
Agarose gel 1.7%TAE

Owi ExactLadder® DNA 100bp is a Ready to Load molecular weight marker . It is designed for the precise sizing of DNA fragments from 100 bp to 1,500 bp.

- **11 sharp bands:** 100 to 1,500 bp
- **Semi-Quantitative:** assess PCR product yield
- **Ready to Load:** load directly onto gel
- **Includes 6X Loading Buffer**
- **Dual tracking dyes:** Orange G and Xylene Cyanol FF
- **Compatible with** Owi EmeralDye Safe DNA Stain

Références	Name	Size
258500	Owi ExactLadder® DNA 100bp Ready to Load	5 lanes (25 µL)
258501	Owi ExactLadder® DNA 100bp Ready to Load	100 lanes (500 µL)