

AS300 Multiple Vacuum Filtration Apparatus

The stainless steel manifold is fitted for three or six filtration units.

Overview

The stainless steel manifold is fitted for three or six filtration units. The apparatus can be autoclaved and sterilized by dry heat at up to 180°C. Suitable only for vacuum operation. If flushing tubes are used, do not exceed 1.3 bar (300 mbar over-pressure). Applications includes: Microbiology quality control, residue analyses, serial filtration carried out rapidly and easily with a common drainage outlet.

Product specifications

AS300/3 Three-place vacuum filtration system, stainless steel 500 ml

Catalog code	10445830

Whatman™ AS310/3 Three-place vacuum filtration system, stainless steel 500 ml, with rapid closure clamp (1/Pk)

Catalog code	10445835
Pack size	1/Pk
Filter Volume	500 ml
Order Information	Recommended for Microbiology Monitors and Alytical Funnels

AS300/5 Three-place vacuum filtration system, stainless steel 100 ml

Catalog code	10445850

DOMINIQUE DUTSCHER SAS