


S·C·A·T®

# The Safest Air Valve for HPLC!






## Best Protection

Solvents can escape, even through the smallest openings. The S.C.A.T. System protects you and your environment from escaping vapour.


## Stable Mixtures

Because the individual components of your solvent mixture can no longer evaporate, your analytic process always remains stable.

## Clean HPLC

Contamination stemming from the surroundings can detrimentally affect your results. The integrated PTFE filter ensures for clean solvents - all the time.

www.scat-europe.com




# Clever!

## Simply turn it around and use the labelling field.

Would you like to note down the expiry date of your valve? No problem - simply turn the clip around and use the labelling field.

### 1 Luer-Lock Adapter
For simple gas injection or for drying tubes.

### 2 PTFE-Filter
Keeps solvents clean.

### 3 Valve Membrane
Blocks solvent vapour.
Integrated check valve prevents the evaporation of solvents.

### 4 Exchange Indicator
Starts working when the button is pressed: never miss an exchange again! Seen sideways on, for improved visibility in the laboratory!

### 5 Labelling Field
Simply turn the clip around and label individually!

### 6 UNF 1/4“ 28G Thread
Fits on to all safety caps.

DOWNLOAD QUE

# Ordering Information



The air valve is available in various packing units (PU), as is best suitable for your HPLC equipment.

| Part No. | Description | PU |
|---|---|---|
| **317 010** | Air Valve for S.C.A.T. Safety Caps | 1 pc. |
| **397 008** | Air Valve for S.C.A.T. Safety Caps | 8 pcs. |
| **397 010** | Air Valve for S.C.A.T. Safety Caps | 10 pcs. |
| **397 050** | Air Valve for S.C.A.T. Safety Caps | 50 pcs. |
| **397 100** | Air Valve for S.C.A.T. Safety Caps | 100 pcs. |

Do you require a specific quotation for your annual demand, or do you have any questions about the product? Contact us! We will gladly assist you further.


**+49 (0) 61 05 / 30 55 86 - 0**

**+49 (0) 61 05 / 30 55 86 - 99**

**info@scat-europe.com**

**www.scat-europe.com**


S·C·A·T®

**The No. 1 System for Solvent Safety.**


# 200157 / 07102021 / 7



DOMINIQUE DU