We make your Lab
a safer Place!

World Premiere:

Safety Waste Cap
for small Neck Threads






© SCAT EUROPE - REPRODUCTION INTERDITE SANS ACCORD PREALABLE DUTSCHER SAS

Safety Waste Cap for small Neck Threads

Discover our most flexible Waste System!



Fits smallest thread sizes

Use empty solvent bottles for HPLC waste

Turn your empty solvent supply bottles into a safe collection container for HPLC waste!

Suitable for MERCK and Wheaton bottles

Including thread sizes GL/S 40 for MERCK and 38-430 for Wheaton bottles.

Easily customizable

- Connects every HPLC system safely to the waste container.
- Blind plugs included for every connection - for a safely closed system.

Workspace Safety

- Connection for exhaust filter - keeps your lab air clean and protects your health.

3 Tubing connectors

- For 5.0 - 11.5 ID tubing
- Including blind plugs for connectors not in use

4 Capillary connectors

- Standard connections for HPLC capillaries with 1.6mm, 2.3mm and 3.2mm OD
- Including all 3 fitting sizes for each connection



Easily expandable

- Put additional HPLC manifolds on top do double the number of connnections!

Threads
- GL 28
- GL 32
- GL 38
- 38-430
- S 40


DOMINIQUE DUTSCHER SAS

Safety Waste Cap for small Neck Thread Bottles
Ordering Information

A 350 028



B 350 032



C 350 038



D 350 039



E 350 040



F 160 532



DOMINIQUE DUTSCHER SAS

Pic.	Part no.	Description	Thread size	Connections for 1.6/2.3/3.2mm OD tubing	Connections for 5.0-11.5 ID tubing	Unit
A	**350 028**	Safety Waste Cap Mini-LISA, GL 28	GL 45	4x	3x	1
B	**350 032**	Safety Waste Cap Mini-LISA, GL 32	GL 45	4x	3x	1
C	**350 038**	Safety Waste Cap Mini-LISA, 38-430	GL 45	4x	3x	1
D	**350 039**	Safety Waste Cap Mini-LISA, GL 38	GL 45	4x	3x	1
E	**350 040**	Safety Waste Cap Mini-LISA, GL/S 40	GL 45	4x	3x	1
F	**160 532**	HPLC Manifold PTFE, GL 14	GL 14	4x	3x	

Safety Waste Cap for small Neck Thread Bottles
Consumables & Accessories

A 107 063



B 107 059



C 107 061



D 160 502



E 117 808



F 160 506



G 107 620



H 490 335



I 490 336



J 407 986



Pic.	Part no.	Description	Diameter	Thread size	Material	Color	Unit
A	**107 063**	PFA fitting with integrated ferrule	3.2 mm ID	UNF 1/4" 28G	PFA	Blue	5
B	**107 059**	PFA fitting with integrated ferrule	2.3 mm ID	UNF 1/4" 28G	PFA	Violet	5
C	**107 061**	PFA fitting with integrated ferrule	1.6 mm ID	UNF 1/4" 28G	PFA	Green	5
D	**160 502**	Blind plug, for capillary connector	-	UNF 1/4" 28G	PFA	Natural	5
E	**117 808**	Tube connector, curved	5.0-11.5 mm OD	NPT 1/8"	PP	Blue	1
F	**160 506**	Blind plug for tube connector	-	NPT 1/8"	PP	Blue	1
G	**107 620**	Blind plug, for exhaust filter connector	-	GL 14	PTFE	White	1
H	**490 335**	Exhaust filter S, service life 3 months each	-	GL 14	PP	Red	4
I	**490 336**	Exhaust filter M, service life 6 months each	-	GL 14	PP	Red	2
J	**407 986**	Exhaust filter L, service life 12 months each	-	GL 14	PP	Red	1

DOMINIQUE DUTSCHER SAS

SCAT - Europe GmbH

Waldecker Straße 7
64546 Mörfelden-Walldorf
Deutschland

+49 (0) 61 05 / 30 55 86 - 0
+49 (0) 61 05 / 30 55 86 - 99
info@scat-europe.com
www.scat-europe.com

SCAT Lab Safety

1704 East Blvd., Suite 101
Charlotte, NC 28203
USA

+1-800-454-1326
info@scatlabsafety.com
www.scatlabsafety.com

DOMINIQUE DUTSCHER SAS