HIRSCHMANN®

Datenblatt
Kapillarröhrchen

Hersteller

HIRSCHMANN LABORGERÄTE

HAUPTSTRASSE 7-15 • 74246 EBERSTADT

Artikelnummer

9201010	9201310	9201510	9201710	9202310
9201012	9201312	9201512	9201712	9202312
9201015	9201315	9201515	9201715	9202315
9201070	9201360	9201570	9201770	9202370
9201075	9201370	9201575	9201775	9202375
9201080	9201375	9201580	9201780	9202380
9201090	9201380	9201590	9201790	9202390
	9201390			

Artikelbeschreibung

KAPILLARRÖHRCHEN BEIDSEITIG OFFEN

Produktbild



Produkteigenschaften

- Beidseitig offen
- Nicht beschichtet

Maße

Artikelnummer	Länge (mm)
920**60	60
920**70	70
920**75	75
920**80	80
920**90	90
920**10	100
920**12	120
920**15	150

Artikelnummer	Außen ⌀ ca. mm	Innen ⌀ ca. mm
92010**	1,00	0,80
92013**	1,35	0,95
92015**	1,55	1,15
92017**	1,75	1,30
92023**	2,35	1,55

Verpackungseinheiten

- Endung *60, *70, *75, *80, *90 und *10 werden zu 10 x 100 Stück verpackt.
- Endung *12 und *15 werden zu 4 x 250 Stück verpackt.
- 92023** werden alle zu 4 x 250 Stück verpackt.

Anwendungsbereich

Labore, sterile Räume und medizinische Einrichtungen.

Material

Artikel 92010** bestehen aus Borosilikatglas 3.3.

Artikel 92013**, 92015**, 92017**, 92023** bestehen aus Sodaglas.

Materialeigenschaften

	max. Gebrauchstemperatur	Transformationstemperatur
Boro 3.3	500 C°	525 C°
Sodaglas	350 C°	550°C ± 10°C

Geltende Normen und Richtlinien

- ISO 9001:2015

Lagerung

Die Produkte sind geeignet für eine Lagerung unter folgenden Voraussetzungen:

- Trockene Lagerung
- Keine übermäßige Hitze
- Kein Gefrieren
- Schutz vor Wetter

DOMINIQUE DUTSCHER SAS

Germany
Hirschmann Laborgeräte GmbH & Co. KG
Hauptstraße 7-15
74246 Eberstadt
Phone +49 7134 511-0

Für weitere Informationen, besuchen Sie gerne: www.hirschmannlab.de

USA
Hirschmann, Inc.
Laboratory Equipment
2700 Holloway Road, Suite 104
Bluegrass Industrial Park
KY 40299 Louisville
Phone: +1 502 2400287

DOMINIQUE DUTSCHER SAS

Die im obenstehenden Datenblatt dargestellten Informationen repräsentieren Standardwerte. Hirschmann übernimmt keine Verantwortung für deren Richtigkeit. Nutzer sind selbst verantwortlich für die Eignung der Produkte, die Durchführung von Voruntersuchungen und die Einhaltung gesetzlicher Bestimmungen. Eine Verwendung der Produkte für Implantationen oder Injektionen ist nicht gestattet. Technische Daten können ohne vorherige Ankündigung geändert werden, dieses Recht behält Hirschmann sich vor.
©Hirschmann Laborgeräte GmbH & Co. KG

HIRSCHMANN®

Product data sheet

Capillary tubes

Manufacturer

HIRSCHMANN LABORGERÄTE

HAUPTSTRASSE 7-15 • 74246 EBERSTADT

Item number

9201010	9201310	9201510	9201710	9202310
9201012	9201312	9201512	9201712	9202312
9201015	9201315	9201515	9201715	9202315
9201070	9201360	9201570	9201770	9202370
9201075	9201370	9201575	9201775	9202375
9201080	9201375	9201580	9201780	9202380
9201090	9201380	9201590	9201790	9202390
	9201390			

Item description

CAPILLARY TUBES BOTH ENDS OPEN

Product image



Product features

- Both ends open
- Not coated

Dimensions

Item number	length (mm)
920**60	60
920**70	70
920**75	75
920**80	80
920**90	90
920**10	100
920**12	120
920**15	150

Item number	Outer ⌀ approx. mm	Inner ⌀ approx. mm
92010**	1,00	0,80
92013**	1,35	0,95
92015**	1,55	1,15
92017**	1,75	1,30
92023**	2,35	1,55

Packaging units

- Suffixes *60, *70, *75, *80, *90 and *10 are packaged in units of 10 × 100 pieces.
- Suffixes *12 and *15 are packaged in units of 4 × 250 pieces.
- 92023** are packaged in units of 4 × 250 pieces.

Area of application

Laboratories, sterile rooms and medical facilities.

Material

Items 92010** are made of borosilicate glass 3.3.

Items 92013**, 92015**, 92017**, 92023** are made of soda-lime glass.

Material properties

	Max. admissible operating temperature	Transformation temperature
Boro 3.3	500 C°	525 C°
Sodaglas	350 C°	550°C ± 10°C

Applicable standards and regulations

- ISO 9001:2015

Storage

The products are suitable for storage under the following conditions:

- Dry storage
- No excessive heat
- No freezing
- Protection from weather conditions

DOMINIQUE DUTSCHER SAS

Germany
Hirschmann Laborgeräte GmbH & Co. KG
Hauptstraße 7-15
74246 Eberstadt
Phone +49 7134 511-0

For more information, please visit:
www.hirschmannlab.de

USA
Hirschmann, Inc.
Laboratory Equipment
2700 Holloway Road, Suite 104
Bluegrass Industrial Park
KY 40299 Louisville
Phone: +1 502 2400287

DOMINIQUE DUTSCHER SAS

The information presented in the above data sheet represents standard values. Hirschmann assumes no responsibility for their accuracy. Users are solely responsible for the suitability of the products, the conduct of preliminary investigations, and compliance with legal regulations. Use of the products for implantations or injections is not permitted. Technical specifications may be changed without prior notice; Hirschmann reserves this right.
©Hirschmann Laborgeräte GmbH & Co. KG