ClearLine	Technical data sheet Serum	Ref : 509509

Iron supplementation treatment on sera

For iron (minus last digit) *E.g.: S1810 => S181R* at the end of the regular CAT n° serum

Analyses:
Tests and results are recorded on CoA.

Treatments:
A 1% Ferric citrate sterile solution is added to the sterile serum at 0.246% (v/v).
Then the serum is well mixed.
Transferrin is normally 20 to 40% bound to iron.
Iron concentration in treated product must be > 500 µg/100ml.

Remarks:
All other parameters than the ones mentioned above remain equal and can be find in each regular (non-treated) serum TDS.

DOMINIQUE DUTSCHER SAS