DOMINIQUE DUTSCHER SAS



















DOMINIQUE DUTSCHER SAS









# CATALOG




2020
new factory running



| | | | | |
|---|---|---|---|---|
| |  |  |  |  |
| **Part No.** | P-1.0-SQ-96 | P-1.6-SQ-96 | P-1.6-SQV-96 | P-2.2-SQ-96 |
| **Description** | 96 square well | 96 square well | 96 square well | 96 square well |
| **Well** | 96 | 96 | 96 | 96 |
| **Volume** | 1.0mL | 1.6ml | 1.6ml | 2.2mL |
| **Well Shape** |  |  |  |  |
| **Size** | 127.8*85.5*24.4mm | 127.8*85.5*31.6mm | 127.8*85.5*31.6mm | 127.8*85.5*41.6mm |
| **Packing size** | 100pcs/case | 96pcs/case | 50pcs/case | 96pcs/case |

| | | | | |
|---|---|---|---|---|
| |  |  |  |  |
| **Part No.** | P-0.5-SQV-96 | P-2.2-SQV-96 | P-2.2-SQVT-96 | P-10-SQV-24 |
| **Description** | 96 square well | 96 square well | 96 square well | 24 square well |
| **Well** | 96 | 96 | 96 | 24 |
| **Volume** | 0.5mL | 2.2mL | 2.2mL | 10ML |
| **Well Shape** |  |  |  |  |
| **Size** | 127.5*85.35*15.1mm | 127.8*85.5*44.2mm | 127.8*85.5*44.2mm | 127.8*85.5*44.2mm |
| **Packing size** | 100pcs/case | 96pcs/case | 96pcs/case | 50pcs/case |

| | | | | |
|---|---|---|---|---|
| |  |  |  |  |
| **Part No.** | P-2.2-SQG-96 | P-2.3-SQV-96 | P-4.6-SQ-48 | P-15-SQV-24 |
| **Description** | 96 square well | 96 square well | 48 square well | 24 square well |
| **Well** | 96 | 96 | 48 | 24 |
| **Volume** | 2.2mL | 2.3ml | 4.6mL | 15ml |
| **Well Shape** |  |  |  |  |
| **Size** | 127.8*85.5*44.2mm | 127.8*85.5*41.6mm | 127.8*85.5*44.2mm | 127.8*85.5*63.2mm |
| **Packing size** | 96pcs/case | 50pcs/case | 96pcs/case | 50pcs/case |

| | | | | |
|---|---|---|---|---|
| |  |  |  |  |
| **Part No.** | P-MT-8U | P-MT-24L/V | P-MT-96V | P-C |
| **Description** | 8 tip comb | 24 tip comb | 96 tip comb | cap |
| **Well** | 8 | 24 | 96 | 0 |
| **Bottom** | U | V | V | / |
| **Assorted** | P-2.2-SQG-96<br>P-2.2-SQ-96 | P-10-SQV-24<br>P-15-SQV-24 | P-2.2-SQV-96<br>P-2.2-SQVT-96 | all deep well plate |
| **Material** | PP | PP | PP | PS |
| **Packing size** | 200pcs/case | 50pcs/case | 50pcs/case | 200pcs/case |

DOMINIQUE DUTSCHER SAS