Lenz Laborglas
GmbH & Co. KG
Am Ried 8
97877 Wertheim
Telefon: (0 93 42) 96 09 0
Fax: (0 93 42) 96 09-30
info@lenz-laborglas.de
www.lenz-laborglas.de

Lenz Laborglas GmbH & Co. KG | Am Ried 8 | 97877 Wertheim

To whom it may concern

Ihre Nachricht/Ihr Zeichen	Unser Zeichen	Name/Durchwahl	Datum
		Harald Karl	2026

DOMINIQUE DUTSCHER SAS

Quality Certificate

We herewith certify, that the glassware manufactured by us are made of Duran borosilicate glass 3.3 and all technical, physical, chemical and thermal properties applying as indicated in the relevant catalogue.

We produce according to DIN 12242/ISO 383, each production step is controlled by a quality management system, based on DIN ISO 9001/EN 29001.

Production on microprocessor controlled half and fully-automatic glass processing machines, particularly supported by robotics and a staff of highly skilled workers, guarantee a continuous high quality standard.

Lenz Laborglas GmbH & Co. KG

Lenz
Laborglas GmbH & Co KG
Am Ried 8
97877 Wertheim/Main

Harald Karl
-Sales Director-