


DOMINIQUE DUTSCHER

# Midi Plus

## Ergonomic Large-Volume Pipetting

Simplifying Progress

SARTORIUS


# Easy and Ergonomic Pipetting of Large Volumes from Bottles or Tubes

**Ergonomic design protects you from RSI**
Midi Plus is a light-weight electronic pipette controller, which allows easy aspiration from bottles and tubes. Its ergonomic design helps to minimize arm and hand raises that are often required while using serological or volumetric pipettes. In addition the electronic operation of Midi Plus reduces the needed finger force compared to a pipette filler with a squeezable bulb. These qualities ease the workload in long periods of pipetting and protect you from Repetitive Strain Injury (RSI).

**Compatible with a wide range of pipettes and tips**
- The Midi Plus is compatible with commonly used glass or plastic pipettes in the 1 – 100 ml volume range.
- It is also compatible with Sartorius 5 ml Optifit Tips with adapter set, available separately.

**Hydrophobic filters prevent over-aspiration**
- The hydrophobic Sartorius filter prevents over-aspiration, which saves the pipetting controller from contamination.
- Filters are available as sterile and non-sterile, to match your purity requirements.

**Save bench space by using the fold-out stand**
- The integrated fold-out stand saves bench space.
- While Midi Plus is resting on the stand, the glass or plastic pipette is kept detached from the table surface, which prevents contamination.




**Find out more**
For more information, please visit
www.sartorius.com

DOMINIQUE DUTSCHER SAS




### Stepless speed control for accurate dispensing

- Aspiration and dispensing speeds can easily be controlled by using the adjustment wheel or by adjusting finger pressure on the operating buttons.
- The speed can be selected according to your application and pipette volume.
- Possibility to dispense by the force of gravity for slow, drop-by-drop dispensing. Just by applying a very small force on the dispensing button.

## Ordering Information and Accessories

| Order Code | Item | Qty \| Unit |
|---|---|---|
| 710931 | Midi Plus | 1 |
| LH-7129120 | Replacement filter, 0.45 µm, non-sterile, autoclavable | 5 |
| LH-7129131 | Replacement filter, 0.2 µm, pre-sterilized | 5 |
| LH-711019 | Adapter set (standard) including nose-cone and the silicone adapter | 1 |
| LH-711017 | Adapter set for 5 ml pipette tip, autoclavable, including nose-cone and the silicone adapter | 1 |
| 780300 | Optifit Tip 5 ml | 100 |
| 780308 | Optifit Tip 5ml | 1000 |

DOMINIQUE DUTSCHER SAS



**Germany**
Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen
Phone +49 551 308 0

**USA**
Sartorius Corporation
3874 Research Park Drive
Ann Arbor, MI 48108
Phone +1 734 769 1600

For further contacts, visit
www.sartorius.com

DOMINIQUE DUTSCHER SAS

Specifications subject to change without notice.
Copyright Sartorius Lab Instruments GmbH & Co. KG.
Status: 09 | 2026 Midi-Plus-Brochure-en-Sartorius