BARRIER®





- Shoulder pleats allow freedom of movement for the upper body



- Forward facing sleeves correspond with natural working position to allow comfort and freedom of movement




- Long sleeve increases freedom of movement and reduces risk of glove sliding

# BARRIER® surgical gown – Fluid Protection Plus


- Lightweight, breathable material for thermal comfort
- Micropores combat perspiration, while polyethylene film prevents fluid penetration
- Impermeable to fluids over the entire surface area
- Ergonomic fit, for greater freedom of movement
- Suitable for procedures with high fluid levels
- Available in High Performance Gowns[1]

DOMINIQUE DUTSCHER SAS

MÖLNLYCKE HEALTH CARE

| High Performance | Ref.No | Description | Size | Length | Width | Disp/Trp |
|---|---|---|---|---|---|---|
|  | 670102 | SG HP, two towels, wrapped | L | 121 cm | 154 cm | 32/64 |
| | 670103 | SG HP, two towels, wrapped | L-L | 131 cm | 154 cm | 30/60 |
| | 670104 | SG HP, two towels, wrapped | XL | 132 cm | 174 cm | 28/56 |
| | 670105 | SG HP, two towels, wrapped | XL-L | 145 cm | 174 cm | 28/56 |
| | 670107 | SG HP, two towels, wrapped | XXL-L | 148 cm | 184 cm | 26/52 |
| | 670106 | SG HP, two towels, wrapped | XXL-XL | 161 cm | 184 cm | 28/56 |

**SG** = Surgical Gown **HP** = High Performance




Measurements for size L

| | |
|---|---|
| **Gown materials** | Gown Body & Sleeve: Polyethylene film and polyester/polyethylene nonwoven. |
| **Product classification** | High performance gowns: European Standard EN 13795:2011 High Performance and ANSI/AAMI PB70:2012 level 3. [1,2] |
| **Packaging** | Single packed gowns provided in a double box packaging system, ensures protection against dust, moisture and damage during transport and storage. Vacuum packaging optimises handling and storage space. |
| **Environment** | At Mölnlycke Health Care we have a long-standing commitment to make our activities have as little effect as possible on the environment, without compromising the efficiency and reliability of our products and services. A sound environmental practice is a priority in our production and at our sites around the world. |

References
1. European Standard EN 13795:2011 Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment – General requirements for manufacturers, processors and products, test methods, performance requirements and performance levels. 2. ANSI/AAMI PB70:2012 Liquid performance and classification of protective apparel and drapes for use in health care facilities.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00.

The Mölnlycke Health Care and BARRIER® trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2015. Mölnlycke Health Care AB. All rights reserved.HQSU0061

MÖLNLYCKE HEALTH CARE

DOMINIQUE DUTSCHER SAS