1903 Protective Gown

Isolation Gown

Product details

Size : L
Colour : Green
Sterile : Non-sterile
Brand : BARRIER

Images



Delivered items

1903-20

Sales released in: Algeria, Austria, Bahrain, Belgium, Bulgaria, Congo (the Democratic Republic of the), Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hong Kong, Hungary, Ireland, Italy, Japan, Latvia, Lebanon, Lithuania, Luxembourg, Malaysia, Montenegro, Netherlands, New Zealand, Norway, Poland, Portugal, Qatar, Romania, Russian Federation, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland
Country of origin: Cambodia
Shelf life: 5 years
Sterilization method: Non-sterile
Production Responsibility: Copious (Cambodia) International Inc., National Road 2, Kilometer 68, Rung Village, Lum Chong Commune, Somrorng District, Takeo Province, Cambodia
Packing information: First packaging layer is a plastic liner/bag. Second layer is a corrugated board transport box.
Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Transport box	50	7323190049721	304x330x330 mm	33.1 dm3	3.4 / 2.9 kg
Pallet	1500	7323190049714	800x1200x1675 mm		

1903-31

Sales released in: Algeria, Austria, Bahrain, Belgium, Bulgaria, Congo, Cyprus, Czechia, Denmark, Estonia, Finland,

DOMINIQUE DUTSCHER SAS

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

France, Germany, Greece, Hong Kong, Hungary, Ireland, Italy, Japan, Latvia, Lebanon, Luxembourg, Mexico, Montenegro, Netherlands, New Zealand, Norway, Poland, Portugal, Qatar, Romania, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: Cambodia

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a plastic liner/bag. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430322067			
Transport box	50	7323190306428	304x330x250 mm	25.1 dm3	2.4 / 2.0 kg
Pallet	2100	7313661085730	800x1200x1894 mm	1818.2 dm3	125.0 kg / -

Material

Animal tissues : No
Human blood derivatives : No
Natural rubber latex : No
Medicinal substances : No
Phthalates : No
Polyvinyl chloride : No

Product Composition Wearing Apparel

Product Component	Composition
Main material	Variant -31: Polypropylene SMMS nonwoven 18 g/m², Variant -20: Polypropylene Spunbond nonwoven 28 g/m²
Neck band	Variant -31: Polypropylene SMMS nonwoven 18 g/m², Variant -20: Polypropylene Spunbond nonwoven 28 g/m²
Elastic	Polyurethane, Polyamide (Nylon)
Tie bands	Variant -31: Polypropylene SMMS nonwoven 18 g/m², Variant -20: Polypropylene Spunbond nonwoven 28 g/m²
Sewing thread	Polyurethane, Polyamide (Nylon)

Technical

Dimension

DOMINIQUE DUTSCHER SAS

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

Dimension reference	Dimension text	Dimension value
A	Length, highest point at shoulder	114 cm
B	Center front	103 cm
C	Chest width	65 cm
D	Top sleeve	56 cm
E	Armhole width	28 cm

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classifications

Regulation type(s)	MDR Class I ns	MDR Class I ns
Intended Purpose :	The intended use is to reduce contamination in a non-sterile environment within health care.	The intended use is to reduce contamination in a non-sterile environment within health care.
MDR Classification Rule :	1	
Conformity Annexes :	IV	
Measuring Function :	No	
Notified body medical devices/PPE :	MPA	MPA
Regulatory Released :	Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Latvia, Lithuania, Luxembourg, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Latvia, Lithuania, Luxembourg, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland

DOMINIQUE DUTSCHER SAS

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001

GMDN Code (Global Medical Device Nomenclature)

35492 Isolation gown, single-use

UNSPSC

42131612 Medical staff isolation or cover gowns

Commodity Code

6210109200 Wearing apparel mainly made of nonwovens

CE Responsibility

Molnlycke Health Care AB, Gamlestadsvagen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000096KF

DOMINIQUE DUTSCHER SAS

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.