Biosigma S.p.A. a Dominique Dutscher Company
Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy |
Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 | E-mail *info@biosigmaeu.com* |
http://www.biosigma.it | http://www.biosigma.com |

SCHEDA TECNICA PRODOTTO / *TECHNICAL DATA SHEET*

Tubi per Drosophila
Drosophila Vials





Descrizione:
Tubi specificamente progettati per il mantenimento di *Drosophila melanogaster*. Realizzati in materiale plastico di alta qualità, garantiscono un'ottima visibilità degli insetti e del substrato di coltura. I tubi sono compatibili con tappi in cotone, spugna o plastica. Ideali per un uso di laboratorio.

Applicazioni:
- Allevamento e mantenimento di ceppi di *Drosophila melanogaster*
- Studi di genetica, biologia dello sviluppo e biologia molecolare
- Esperimenti di incrocio e osservazione fenotipica
- Attività di ricerca e didattica in laboratori accademici e industriali

Materiali:
- PS trasparente (non sterilizzabile),
- PP traslucido (sterilizzabile)

Caratteristiche:
Disponibili in due formati:
- Stretto (25x95 mm)
- Largo (28,5x95mm)

Confezione:
Vassoi da 100 (10x10), in confezione da 500 pz, o sfusi

Tappi:
Disponibili anche tappi in TPV da 25 mm o 28,5 mm, confezionati in sacchetti da 10 x 100 pezzi

Description:
Vials designed specifically for the breeding, and storage of Drosophila melanogaster. Made of high-quality plastic material, they ensure excellent visibility of both the insects and the culture medium. The vials are compatible with cotton, sponge, or plastic plugs. Ideal for laboratory use.

Applications
- *Breeding and maintenance of Drosophila melanogaster strains*
- *Genetic, developmental biology, and molecular biology studies*
- *Crossing experiments and phenotypic observation*
- *Research and teaching activities in academic and industrial laboratories*

Materials
- *Transparent PS (non-sterilizable)*
- *Translucent PP (sterilizable)*

Features
Available in two formats:
- *Narrow (25 × 95 mm)*
- *Wide (28.5 × 95 mm)*

Packaging
100-unit trays (10 × 10), packed in boxes of 500 pcs, or bulk packed.

Caps:
TPV caps are also available in 25 mm and 28.5 mm sizes, packed in 10 × 100 pcs bags

SCHEDA TECNICA PRODOTTO / *TECHNICAL DATA SHEET*

INFORMAZIONI PER L'ORDINE / *ORDERING INFORMATION*

Codice *Code*	Materiale *Material*	Ø x h (mm)	Versione / *Version*	Unità vendita (pz) *Sales unit (pcs)*
789001B	PS	25 x 95	Vassoio da 100 pz *100 pc tray*	500
789002	PS	28.5 x 95	Vassoio da 100 pz *100 pc tray*	500
789001PP	PP	25 x 95	Vassoio da 100 pz *100 pc tray*	500
789002PP	PP	28.5 x 95	Vassoio da 100 pz *100 pc tray*	500
789008	PS	25 x 95	Sfusi / *Bulk*	500
789009	PS	28.5 x 95	Sfusi / *Bulk*	500
789008PP	PP	25 x 95	Sfusi / *Bulk*	500
789009PP	PP	28.5 x 95	Sfusi / *Bulk*	500

*TAPPI/**CAPS***

Codice *Code*	Diametro *Diameter*	Materiale	Unità vendita (pz) *Sales unit (pcs)*
BFF025	25 mm	Gomma/*Rubber*	10 x 100
BFF028	28,5 mm	Gomma/*Rubber*	10 x 100