AUXILAB
Material de laboratorio
Laboratory supplies

PRODUCT COMPLIANCE STATEMENT

Description: KNOB LID DESICCATOR, 200 mm
Reference: BDD003
Brand: Endo Glassware

Technical information:

Inner diameter	Outer diameter	Lid height	Total height
200 mm	270 mm	105 mm	280 mm

Material composition:

The product is manufactured from glass with the following chemical composition:

SiO_2 / Na_2O / CaO

Compliance statement:

AUXILAB S.L. hereby certifies that the product identified above complies with the applicable product specifications, based on the manufacturer's documentation.

The product:

- Contains not intentionally added PFAS and complies with applicable PFAS restrictions
- Has been manufactured under a Quality Management System compliant with ISO 9001

Issued by: Auxilab
Place: Beriáin (Navarra), Spain
Date: September 23th, 2026

AUXILAB, S. L.
N.I.F. ESB31072218
MATERIAL PARA LABORATORIO
Polig. Morea Norte, 6 - Tel. 34-948 310 513
31191 BERIAIN - NAVARRA (ESPAÑA)