PP Heavy-Duty Vacuum Bottle

The perfect container for your applications

WATCH VIDEO

121°C / 1bar / 20 min Autoclavable

Sturdy and Shatterproof

- Thick-walled construction withstands full vacuum for 24 hours.*
- Made from rugged polypropylene (PP), an excellent alternative to glass.

* Test at room temperature (25°C).

High Quality

- EPDM and silicone o-ring guaranteed leak proof under vacuum conditions.

Volume Identification

- Volume graduations are marked at side of bottle for easy identification.

Chemical Resistant & Autoclavable

- Good resistance to a wide range of chemicals.
- Autoclavable at 1 bar and 121°C for 20 minutes.

BPA FREE

FDA APPROVED



Applications

- Vacuum applications
- Waste bottle, safety bottle
- Aseptic liquid transfer bottle
- HPLC solvent container

DOMINIQUE DUTSCHER SAS