Codes: 200320, 200321, 200322
Autoclave resistants red bags (135ºC)

Catalogue page 2007: 97

Code	Dimensions		Volume		Case quantity	Case weight (kg)	Case volume (m^3)
	centimeters	inches	liters	gallons			
200320	48 x 58	19 x 23	24.37	6.44	200	5.10	0.011
200321	65 x 91	25 x 35	75.10	19.84	200	10.14	0.020
200322	78 x 96	31 x 38	120.98	31.96	200	13.25	0.029



SPECIFIC USE AND DESCRIPTION

Made of high molecular weight polypropylene 0.055 mm thick.

Bags printed with biohazard warning symbol and precautionary procedures in four languages: English, Spanish, French and German. They feature a sterilization indicator patch that darkens when exposed to steam sterilization.

DOMINIO DEUTSCHER S.A.S

Elaborated	Verified	Approved

INSTRUCTIONS FOR USE IN AUTOCLAVE

Autoclavables up to **135ºC** (275ºF), during **20 minutes** at sustained heat level. At this temperature, autoclave saves time as the standard process at 121º C requires 30 minutes to ensure sterilization.

Attention:
Do not close bags hermetically in the autoclave process.

STANDARDS

ASTM D 1709-01 *Dart-Impact-1000 grams.*

ASTM D 1922-00a *Elmendorf tear Resistance*
- MD – 200 grams
- TD – 1,300 grams

DOMINIQUE DUTSCHER SAS

Elaborated	Verified	Approved