**Biosigma S.p.A.**
**to Dominique Dutscher Company**
Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |
Fax ++39 0426 302228 | SMS ++39 348 4077376 |
E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

Data Update
27/09/2022

**SCHEDA TECNICA PRODOTTO / *TECHNICAL DATA SHEET***

# Eutectic plates




## Cooling elements for temperature-controlled transport

| Specific | Code | BS4008 | BS4009 | BS4010 |
|---|---|---|---|---|
| | Dimensions | 160x110x14 mm | 250x140x14 mm | 320x240x14 mm |
| | Weight | 200 g | 410 g | 820 g |

| | | | |
|---|---|---|---|
| **Description** | Ultra-flat eutectic plates for the transport of biological samples, drugs and blood components | | |
| **Features** | • Materials | •Slab | Polyethylene |
| | | •Liquid | Mixture of water and salts |
| | | All components are latex free | |
| | • Plates should be frozen -15°C/-18°C for at least 24 hours before use | | |

| | | |
|---|---|---|
| **Cleanliness disinfection and sterilization** | **•Cleanliness** | Neutral surfactant-based detergents |
| | **•Disinfection** | pH neutral disinfectants |
| | **•Sterilization** | NOT autoclavable |
| | **ATTENTION** | For cleaning and disinfection, absolutely avoid the use of:<br>• Acetone • Solvents • Thinners • Alcohol<br>• Concentrated acid solutions<br>• Concentrated basic solutions<br>• Sodium hypochlorite solutions |
| **Maintenance** | Every 300 cycles of use or once a year, inspect the product for integrity.<br>In case of wear, replace the plate. | |
| **Guarantee** | • The product is covered by a 24-month warranty from the date of delivery to the customer.<br>• During the warranty period, the free replacement of the product is ensured if the malfunction/defect is attributable to production defects or material anomalies.<br>• Free replacement is not guaranteed in the event that the malfunction/defect is attributable to improper use of the product or failure to comply with the instructions for use. | |

DOMINIQUE DUTSCHER SAS