BAGFILTER
PIPET & **ROLL**

Filter bags with pipet access & built-in closure



The smart sticker bag!

DOMINIQUE DUTSCHER SAS

interscience

Great lab work efficiency!

PATENTED

BagFilter® Pipet & Roll bags, with **BagRack® *Slide*** storage rack, are an innovative solution for ergonomic work. They are most handy for pathogens research, such as *Listeria* and *Salmonella* analyses, when repeat pipetting from the sample bag is required.

1 Roll

Simply **roll** your bag **3 times** and **close it** with the pre-cut **sticker** for a **waterproof** extra-resistant closure.

2 Slide

Easily **slide your bag** from the side to store your bags during multiple pipettings. Enjoy great **ergonomics** and a **safe**, **cross-contamination free** work environment.




3 Pipet

Open/close the pipet access easily with the **reclosable sticker**. The filtrate is debris-free thanks to the filter. When multiple pipettings are required for **Pathogens research**, access with pipet tips is made easy by the **lower pipetting height**, without risks of musculoskeletal disorders.

DOMINIQUE DUTSCHER SAS

Quick sticker closure



Without any extra tools, simply close with the built-in closure sticker.

BAGFILTER® **PIPET** & **ROLL**


Filter bag with built-in closure

With other products

Bag + Closing system


Extra resistant closure

Get an extra-resistant closure within seconds. Waterproof and resistant up to 5 kg pressure.



BAGFILTER® **PIPET** & **ROLL**	With other products
5 kg	2 kg

Easy pipetting access



During multiple pipettings for pathogens research, open and close the pipet access easily without opening the bag from the top.

Number of manipulations


BAGFILTER® **PIPET** & **ROLL**

1

With other products

1 > 2 > 3 > 4


Slide your bag



Use the BagRack® *Slide* as a stand for your multiple pipetting and sample storage, for great ergonomic work. Prevents R.S.I. (Repetitive Strain Injury)

BAGFILTER® **PIPET** & **ROLL**	With other products
12 s	24 s


Pathogens research



The low pipetting height is great for work under laminar flow cabinets and ideal with pipet tips for Pathogens research.

BAGFILTER® **PIPET** & **ROLL**	With other products
10 cm	30 cm

Rapid & hygienic disposal



After work, waste disposal is much quicker and without any sample dripping.

Quick

- BAGFILTER® **PIPET** & **ROLL**
- With other products

Clean

- BAGFILTER® **PIPET** & **ROLL**
- With other products

DOMINIQUE DUTSCHER SAS

3 models for your applications



BAGFILTER®
PIPET & ROLL



BAGFILTER®
PIPET



BAGFILTER®
ROLL

Built-in closure sticker	✓	-	✓
Reclosable easy pipetting sticker	✓	✓	-

Reference	**111 720**	**111 710**	**111 700**
Type of bag	Bag with lateral filter		
Material of the bag	MultiLayer®: reinforced multicoated complex		
Type of filter	Non-woven		
Filter porosity	< 250 microns		
Max volume	400 mL		
Optimal blending volume	50 - 300 mL		
Optimal pipetting volume	250 mL		
Bag dimensions	190 x 300 mm		
Gamma ray treatment	5 to 12 KGy		
Shelf-life	5 years		
Resistance temperature	From -20°C to 45°C (48h) / 80 °C (4h) - From -4°F to 113°F (48h) / 176°F (4h)		
In compliance with	ISO 7218 - ISO 6887- FDA-BAM standards		
Pack of	25 bags		
Box of	500 bags		

Ask for free samples

Certified manufacturing
Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.

MADE IN FRANCE

Complementary products



Dilu***Flow***®
Gravimetric dilutors
Ref 503 205



BagMixer®
Lab blenders
Ref 025 100



BagRack® *Slide*
Storage rack for bags
Ref 221 440

Your local distributor

www.interscience.com

interscience

PARIS
Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770 - Email: sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

SHANGHAI
Phone: +86 (0)21-64739390 - Email: sales.china@interscience.com

SINGAPORE
Phone: +65 6977 7232 - Email: sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

5592 BagFilterPipet&RolL B 0322. Pictures and information in this brochure are not contractually binding. For presentation reasons, the pipetting zone is shown on the left of the drawing. BagFilter®, DiluFlow® and BagMixer® are trademarks. This document is the translation of the french brochure and is not contractually binding. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860. Printed by Accent Tonic.

EN

DOMINIQUE DUTSCHER SAS