eppendorf

NEW accessory for all Eppendorf pipettes



Make Your Mark

Label your pipettes with ColorTag marking rings

Make your pipette your own and keep track of essential information with our new ColorTag silicone pipette marking rings, available in six sizes and twelve colors. ColorTag marking rings attach easily to any Eppendorf pipette and many other liquid handling instruments, offering color-coded identification from afar. Customize them with pens* for personalized labeling!

Features

> Flexibly label any pipette for personal use or specific applications.
> Prevent errors from pipette misidentification.
> Add color and joy to your lab environment!
> Autoclavable and extensively tested for lab use.


*Tested with ballpoint pens, felt tip and/or permanent markers. Soak/wipe with ethanol to remove ink (not suitable when using permanent markers).

Learn more
www.eppendorf.link/ColorTag

DOMINIQUE DUTSCHER SAS

Ordering information

ColorTag, pipette marking rings for Eppendorf and third party liquid handling instruments

Color	Eppendorf fit (fits pipette lower part, unless specified otherwise)	1-channel pipettes up to 1 mL	Upper part of all pipettes > 2 mL pipettes	> 5 mL pipettes > Easypet® 3, aspirating cone	> 10 mL pipettes > All Multipette® dispensers > Easypet® 3, grip > Pipette Carousel 2 / Pipette Stand 2	> 8- or 16-channel pipettes > Eppendorf Top Buret > Varispenser® 2 / 2x	> 12- or 24-channel pipettes > Move It® pipettes
	Size (Ø in mm)	19	24	27	34	50	73
	Pieces per pack	10	10	5	5	5	5
light blue		3102 660 000	3102 661 000	3102 662 000	3102 663 000	3102 664 000	3102 665 000
light green		3102 660 001	3102 661 001	3102 662 001	3102 663 001	3102 664 001	3102 665 001
light yellow		3102 660 002	3102 661 002	3102 662 002	3102 663 002	3102 664 002	3102 665 002
light orange		3102 660 003	3102 661 003	3102 662 003	3102 663 003	3102 664 003	3102 665 003
light pink		3102 660 004	3102 661 004	3102 662 004	3102 663 004	3102 664 004	3102 665 004
light violet		3102 660 005	3102 661 005	3102 662 005	3102 663 005	3102 664 005	3102 665 005
neon blue		3102 660 010	3102 661 010	3102 662 010	3102 663 010	3102 664 010	3102 665 010
neon green		3102 660 011	3102 661 011	3102 662 011	3102 663 011	3102 664 011	3102 665 011
neon yellow		3102 660 012	3102 661 012	3102 662 012	3102 663 012	3102 664 012	3102 665 012
neon orange		3102 660 013	3102 661 013	3102 662 013	3102 663 013	3102 664 013	3102 665 013
neon pink		3102 660 014	3102 661 014	3102 662 014	3102 663 014	3102 664 014	3102 665 014
neon magenta		3102 660 015	3102 661 015	3102 662 015	3102 663 015	3102 664 015	3102 665 015

ColorTag, pipette marking rings sample bag	Pcs. per pack	Order no.
Each ColorTag pipette marking ring size in two colors (all 12 colors included)	12	3102 666 000

Your local distributor: **www.eppendorf.link/contact**
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, Easypet®, Multipette®, Varispenser®, and Move It® are registered trademarks of Eppendorf SE, Germany. All right reserved, including graphics and images.
Copyright © 2025 by Eppendorf SE. Order No: AA01 049 320 – EN1/0325/BE

DOMINIQUE DUTSCHER SAS