S.C.A.T. Europe Sets Overview
HPLC Supply and Waste
Generation V2.0

THE System for Solvent Safety

HPLC Starter Kits: Supply

with 4x SafetyCaps incl.:
- 4x Air valves
- Fittings
- Blind plugs



HPLC Supply-Sets: Supply

comprising:
- SafetyCap incl.: Air valve, Fittings, Blind plugs
- Laboratory glass bottle
- Capillary / Capillaries
- Suction filters



HPLC Waste-Sets: Waste

comprising:
- SafetyWasteCap incl.: Fittings, Tube connector, Blind plugs
- Exhaust filter
- Canister



HPLC Supply and Waste Set: Supply and Waste - Complete Set - 1 Set - suitable for all HPLC systems!

comprising:
- 4x SafetyCaps incl. Air valves, Fittings, Blind plugs
- 4x Laboratory glass bottles
- 1x SafetyWasteCap incl. Fittings, Tube connectors, Blind plugs
- 1x Exhaust filter
- 1x Canister



Part No.		Description
399200		HPLC **Starter** Kit **1** - **Supply**, V2.0: **3x SafetyCaps I** (307019), **1x SafetyCap II** (307909), 1x blind plug, 5x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
399201		HPLC **Starter** Kit **2** - **Supply**, V2.0: **4x SafetyCaps II** (307909), 4x blind plugs, 8x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
399202		HPLC **Starter** Kit **3** - **Supply**, V2.0: **4x SafetyCaps III** (307910), 8x blind plugs, 12x fittings 1.6/2.3/3.2 mm, 4x air valves (317010)
307300		HPLC **Supply**-Set **I**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & clear glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter** + air valve (317010)
307301		HPLC **Supply**-Set **I**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **square & clear glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter** + air valve (317010)
307312		HPLC **Supply**-Set **I**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & brown glass + 1x 1.5 m capillary 3.2 mm OD, 1x PTFE suction filter** + air valve (317010)
307303		HPLC **Supply**-Set **II**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & clear glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters** + air valve (317010)
307305		HPLC **Supply**-Set **II**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **square & clear glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters** + air valve (317010)
307313		HPLC **Supply**-Set **II**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & brown glass + 2x 1.5 m capillaries 3.2 mm OD, 2x PTFE suction filters** + air valve (317010)
307304		HPLC **Supply**-Set **III**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters** + air valve (317010)
307306		HPLC **Supply**-Set **III**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **square & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters** + air valve (317010)
307314		HPLC **Supply**-Set **III**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **round & brown glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters** + air valve (317010)
307347		HPLC **Supply**-Set **III**, V2.0: 1x SafetyCap I GL45, 1x 1000 ml Laboratory bottle **Youtility & clear glass + 3x 1.5 m capillaries 3.2 mm OD, 3x PTFE suction filters** + air valve (317010)
307307		HPLC **Waste**-Set **1**, V2.0: SafetyWasteCap GL45, 2x 2.3/3.2 mm OD, 1x 6.4 - 9 mm ID (307923), 5 liter canister, PE-HD (107951), exhaust filter M, with splash guard and exchange indicator (310535)
307310		HPLC **Waste**-Set **2**, V2.0: SafetyWasteCap GL45, 2x 2.3/3.2 mm OD, 1x 6.4 - 9 mm ID (307923), 10 liter canister, PE-HD (107952), exhaust filter M, with splash guard and exchange indicator (310535)
307328		HPLC **Waste**-Set **3**, V2.0: SafetyWasteCap S50, 2x 3.2 mm OD, 1x 5.0 - 11.5mm ID (108025), 5 liter space saving canister, PP (107998), exhaust filter S, with splash guard and exchange indicator (310534)
307337		HPLC **Supply and Waste Set, V2.0, suitable for all HPLC systems**

DOMINIQUE DUTSCHER SA

#200web / 15122020