3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
1 kb DNA Ladder - for DNA electrophoresis

Product Number: **D0428**

Storage Temperature: -20 ℃

TEST	Specification
Agarose Gel Electrophoresis Band pattern consistent with sizes of 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 4.0, 5.0, 6.0, 8.0 and 10.0 kb.	Conforms

Specification: PRD.0.ZQ5.10000069150

DOMINIQUE DUTSCHER SAS

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.