dominique Dutscher

Product code : L0130
Product name : Ham's F10 w/ L-Glutamine w/ 25 mM Hepes
Désignation : Ham's F10 avec L-Glutamine avec 25 mM Hepes

CAS Number	Components	Quantity in g/l
10035-04-8	Calcium Chloride Dihydrate	0.04410000
7487-88-9	Magnesium Sulfate Anhydrous	0.07464000
7758-99-8	Cupric Sulfate Pentahydrate	0.00000250
7782-63-0	Ferrous Sulfate Heptahydrate	0.00083400
7447-40-7	Potassium Chloride	0.28500000
7778-77-0	Potassium Phosphate Monobasic Anhydrous	0.08300000
7647-14-5	Sodium Chloride	6.80000000
7558-79-4	Sodium Phosphate Dibasic Anhydrous	0.15370000
7446-20-0	Zinc Sulfate Heptahydrate	0.00002880
50-99-7	D-Glucose Anhydrous	1.10000000
56-40-6	Glycine	0.00751000
56-41-7	L-Alanine	0.00900000
1119-34-2	L-Arginine Monohydrochloride	0.21100000
5794-13-8	L-Asparagine Monohydrate	0.01501000
56-84-8	L-Aspartic acid	0.01330000
7048-04-6	L-Cysteine Monohydrochloride Monohydrate	0.03500000
56-86-0	L-Glutamic Acid	0.01470000
56-85-9	L-Glutamine	0.14600000
5934-29-2	L-Histidine Monohydrochloride Monohydrate	0.02100000
73-32-5	L-Isoleucine	0.00260000
61-90-5	L-Leucine	0.01310000
657-27-2	L-Lysine Monohydrochloride	0.02930000
63-68-3	L-Methionine	0.00448000
63-91-2	L-Phenylalanine	0.00496000
147-85-3	L-Proline	0.01150000
56-45-1	L-Serine	0.01050000
72-19-5	L-Threonine	0.00357000
73-22-3	L-Tryptophan	0.00060000
69847-45-6	L-Tyrosine Disodium Salt Dihydrate	0.00261000
72-18-4	L-Valine	0.00350000
67-48-1	Choline Chloride	0.00069800
58-85-5	D-Biotin	0.00002400
137-08-6	D-Ca Pantothenate	0.00071500
59-30-3	Folic Acid	0.00132000
87-89-8	Myo-Inositol	0.00054100
98-92-0	Nicotinamide (Nicotinic acid amide)	0.00061500
58-56-0	Pyridoxine Hydrochloride	0.00020600
83-88-5	Riboflavin	0.00037600
67-03-8	Thiamine Hydrochloride	0.00100000
68-19-9	Vitamine B12	0.00136000
7365-45-9	Hepes Free Acid	5.95800000
68-94-0	Hypoxanthine	0.00408000
34487-61-1	Phenol Red Sodium Salt	0.00130000

DOMINIQUE DUTSCHER SAS

dominique Dutscher

Product code : L0130
Product name : Ham's F10 w/ L-Glutamine w/ 25 mM Hepes
Désignation : Ham's F10 avec L-Glutamine avec 25 mM Hepes

CAS Number	Components	Quantity in g/l
113-24-6	Sodium Pyruvate	0.11000000
1077-28-7	Thioctic Acid	0.00021000
50-89-5	Thymidine	0.00073000
144-55-8	Sodium Bicarbonate	1.20000000
WATER		983.61827970

DOMINIQUE DUTSCHER SAS