# Composition

**P04-30500 EBSS**

| | Components | mg/L |
|---|---|---|
| Inorganic Salts | Potassium chloride | 400.00 |
| | Sodium chloride | 6,800.00 |
| | Sodium dihydrogen phosphate x $H_2O$ | 140.00 |
| | Calcium chloride anhydrous | 200.00 |
| | Magnesium sulfate x $7H_2O$ | 200.00 |
| Other Components | D(+)-Glucose anhydrous | 1,000.00 |
| | Phenol red | 10.00 |

DOMINIQUE DUTSCHER SAS