

This drawing and/or any technical information relating to it issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances


A0




ENGRAVING DETAILS

ALL CHARACTERS TO BE RAISED 0.20
NUMBERS TO BE 1.5 HIGH
LETTERS TO BE 2.5 HIGH







SECTION B-B




SECTION C-C

121.75 INT - 4° INC.
125.56 - 4° INC.





M 2:1



VIAL FORM 384 pos.
SCALE 10:1

Code	Description	Plates/Case
4ti-0384	FrameStar® 384well PCR plate, clear wells, purple frame	50
4ti-0384/B	FrameStar® 384well PCR plate, clear wells, blue frame	50
4ti-0384/C	FrameStar® 384well PCR plate, clear wells, clear frame	50
4ti-0384/G	FrameStar® 384well PCR plate, clear wells, green frame	50
4ti-0384/R	FrameStar® 384well PCR plate, clear wells, red frame	50
4ti-0384/X	FrameStar® 384well PCR plate, clear wells, black frame	50
4ti-0385	FrameStar® 384well PCR plate, white wells, black frame	50
4ti-0386	FrameStar® 384well PCR plate, black wells, black frame	50

ISSUE	DATE	AMENDMENT

DO NOT SCALE
IF IN DOUBT ASK

4titude®
www.4ti.co.uk

CUSTOMER:
CUSTOMER DRG No: ISSUE:
MATERIAL: FRAME: Polycarbonate; TUBES: Polypropylene
TOLERANCES: GENERAL TOLERANCES (UNLESS OTHERWISE STATED) 0 ±0.20; 0.0 ±0.10; 0.00 ±0.05; ANGULAR ±0°30' (ALL DIMENSIONS IN MILLIMETRES)

DRAWN: ZOLTAN E.
DATE: 14.10.2014
CHECKED: JAMES B.
DATE:
SCALE:

TITLE: FrameStar 384
Marketing drawing
DRAWING No: 4ti-0384-NF
SHEET No: 1 OF 1 ISSUE: E

DOMINIQUE DUTSCHER SAS