CONDENSYN
WATERLESS REFLUX CONDENSERS


CondenSyn and CondenSyn MINI - Waterless reflux condensers:
- Suitable for reflux and distillation
- Chemical-resistant and durable
- Long-term use up to 230°C
- Short-term use up to 400°C
- No water, no risk of flooding
- Easy to use, no hoses or water flow to control
- Can be used in vacuum / inert conditions
- Suitable for air-sensitive work when connecting
- Lightweight - less clamping required, ideal for angled ports
- Space efficient - no need to be near a water source
- Easy to clean - rinse with solvent/bath solvent
- Visibility - Can show condensation during use and contamination during cleaning
- Can be dried in an oven at standard temperatures (~120°C) or heat gun

NEW 2023 !

Designation	Asynt part no.	Milian part no.	Units/ box	CHF/ box
CondenSyn MINI 200 mm				
Asynt CondenSyn MINI: length 200 mm with B14 connections	**GB-C-200-B14**	**041-412121**	1	**NC -**
Asynt CondenSyn MINI: length 200 mm with B19 connections	**GB-C-200-B19**	**041-412122**	1	**NC -**
Asynt CondenSyn MINI: length 200 mm with B24 connections	**GB-C-200-B24**	**041-412123**	1	**NC -**
Asynt CondenSyn MINI: length 200 mm with B29 connections	**GB-C-200-B29**	**041-412124**	1	**NC -**
CondenSyn 350 mm				
Asynt CondenSyn 350 : length 350 mm with B14 connections	**GB-C-350-B14**	**041-412125**	1	**NC -**
Asynt CondenSyn 350 : length 350 mm with B19 connections	**GB-C-350-B19**	**041-412126**	1	**NC -**
Asynt CondenSyn 350 : length 350 mm with B24 connections	**GB-C-350-B24**	**041-412128**	1	**NC -**
Asynt CondenSyn 350 : length 350 mm with B29 connections	**GB-C-350-B29**	**041-412134**	1	**NC -**
CondenSyn 450 mm				
Asynt CondenSyn 450 : length 450 mm with B14 connections	**GB-C-450-B14**	**041-412138**	1	**NC -**
Asynt CondenSyn 450 : length 350 mm with B19 connections	**GB-C-450-B19**	**041-412139**	1	**NC -**
Asynt CondenSyn 450 : length 350 mm with B24 connections	**GB-C-450-B24**	**041-412140**	1	**NC -**
Asynt CondenSyn 450 : length 350 mm with B29 connections	**GB-C-450-B29**	**041-412141**	1	**NC -**
Asynt CondenSyn 450 : length 350 mm with B34 connections	**GB-C-450-B34**	**041-412142**	1	**NC -**



CONTACT US FOR MORE INFORMATION ON ASYNT
MANY OTHER OPTIONS ARE AVAILABLE
DEMONSTRATION UNITS ARE AVAILABLE



Chemin du Champ-des-Filles n°4A
Campus Spark
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
www.milian.com

DRYSYN CLASSIC
SINGLE-POSITION HEATING BLOCKS

- Made of anodised aluminium
- Chemical resistant and durable
- Reduces the risk of fire (no spills, no toxic fumes)
- Maintenance-free
- Consistent heating times
- Suitable for all hobs
- Clear visibility of reaction
- Heats up to 500°C

Versatile:
- Base with heat-resistant handles
- Base accepts a standard 1,000 ml round-bottomed vial
- Base fits all standard vials up to 25 ml
- Base fits all standard flasks up to 25 ml
- Adjustable feet for full compatibility with all hot plate stirrers



Designation	Asynt part no.	Milian part no.	Units/ box	CHF/ box
DrySyn Single Position Heating: 25 mL - 1000 mL				
DRYSYN COMPLETE KIT WITH BASE 1000ml WITH HANDLES AND 4 ADAPTERS	**ADS1**	**041-411987**	1	**NC -**
DRYSYN BASE 1000ml WITH RESISTENT HANDLES	**ADS1-B**	**041-411988**	1	**NC -**
Accessories				
DRYSYN ADAPTER 25ml	**ADS1-25**	**041-411989**	1	**NC -**
DRYSYN ADAPTER 50ml	**ADS1-50**	**041-411990**	1	**NC -**
DRYSYN ADAPTER 100ml	**ADS1-100**	**041-411991**	1	**NC -**
DRYSYN ADAPTER 250ml	**ADS1-250**	**041-411992**	1	**NC -**
DRYSYN ADAPTER 500ml	**ADS1-500**	**041-411993**	1	**NC -**
DRYSYN MULTI-E BASE WITH HANDLES	**ADS3-B**	**041-412001**	1	**NC -**
DRYSYN MULTI-E PLIERS 3 POSITIONS	**ADS3-C**	**041-412002**	1	**NC -**
DRYSYN MULTI-E ADAPTER 5ml X3	**ADS3-5**	**041-412003**	1	**NC -**
DRYSYN MULTI-E ADAPTER 10ml X3	**ADS3-10**	**041-412004**	1	**NC -**
DRYSYN MULTI-E ADAPTER 25ml X3	**ADS3-25**	**041-412005**	1	**NC -**
DRYSYN MULTI-E ADAPTER 50ml X3	**ADS3-50**	**041-412006**	1	**NC -**
DRYSYN MULTI-E ADAPTER 100ml X3	**ADS3-100**	**041-412007**	1	**NC -**
DRYSYN MULTI-E ADAPTER VIALS X3	**ADS3-R-(size)**	**041-412008**	1	**NC -**
DRYSYN MULTI-E COMPLETE KIT	**ADS3**	**041-412016**	1	**NC -**
DRYSYN MULTI-E SUPER KIT	**ADS3-A**	**041-412017**	1	**NC -**
DRYSYN PARALLEL SYNTHESIS KIT	**ADS3-S-(size)**	**041-412018**	1	**NC -**

CONTACT US FOR MORE INFORMATION ABOUT ASYNT
MANY OTHER OPTIONS AVAILABLE
DEMONSTRATION UNITS AVAILABLE

Asynt
www.asynt.com



Chemin du Champ-des-Filles n°4A
Campus Spark
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
www.milian.com