milian®
DUTSCHER GROUP

eppendorf


GoGREEN

ep.T.I.P.S.® BIOBASED EPPENDORF

- Improve your laboratory's carbon footprint and move to epT.I.P.S.®BioBased with the new Eppendorf Reload. These pipette tips are made from second-generation recycled raw materials, obtained exclusively from waste and residues from vegetable oil production, such as used edible oil.
- Significantly reduce your amount of laboratory waste - use refills rather than disposable racks
- Trust the promise of sustainability - the production site is certified by ISCC PLUS (International Sustainability & Carbon Certification) for traceability throughout the supply chain.
- Pipette tips certified by ACT (Accountability, Consistency, and Transparency), the independent validation body for assessing the sustainability of products
- **BioPure:** Filter-free, sterile and free of pyrogens, RNAs, DNases, DNA, ATP and PCR inhibitors
- **PCR Clean and Sterile:** With filter, sterile and free of pyrogens, RNAs, DNases, DNA and ATP, or PCR inhibitors
- **SealMax:** Self-sealing filter. An innovative two-phase filter system safely retains all liquids and moisture. Sterile and free of pyrogens, RNAs, DNases, DNA and ATP, and PCR inhibitors







Premium material, certified quality and purity:

- Excellent sample integrity thanks to the absence of lubricants, plasticisers and biocides
- Batch-specific purity certificates, issued by an independent laboratory, can be downloaded from the Eppendorf website.
- Global quality certificates for compliance with Eppendorf standards and general quality requirements
- Available in volumes from 10 µL to 1,250 µL, long in refills for use with epT.I.P.S.® Box 2.0
- Available as epT.I.P.S.® BioBased Biopur®, ep Dualfilter T.I.P.S.® BioBased PCR clean and sterile, ep Dualfilter T.I.P.S.® SealMax® BioBased Biopur®

milian®
DUTSCHER GROUP

Chemin du Champ-des-Filles n°4A
Campus Spark
CH-1228 PLAN-LES-OUATES
Tel. +41 (0) 22 884 16 00
www.milian.com

milian DUTSCHER GROUP

eppendorf

GoGREEN

M = Medium
L = Long

ep.T.I.P.S.® BioBased Reloads

Part number	Designation	Volume	Packaging unit
	Autoclavable reusable boxes		
041-033360	epT.I.P.S.® Box 2.0 empty, epQuality, for epT.I.P.S.® Reloads, (0 trays × 0 tips)	10 - 300 µL	1 reusable autoclavable box
041-033361	epT.I.P.S.® Box 2.0 empty, epQuality, for epT.I.P.S.® Reloads, (0 trays × 0 tips)	1000 - 1250 µL	1 reusable autoclavable box
041-033362	epT.I.P.S.® Box 2.0 empty, epQuality, for epT.I.P.S.® Reloads, (0 trays × 0 tips)	1250 µL long - 5 mL	1 reusable autoclavable box
	BioBased tips - BioPur Reloads		
041-033354	epT.I.P.S.® Reloads, BioBased, Biopur, 40 mm, medium grey	0.1 to 20 µL	5 trays of 96 tips (480 tips)
041-033355	epT.I.P.S.® Reloads, BioBased, Biopur, 53 mm, yellow	2 to 200 µL	5 trays of 96 tips (480 tips)
041-033356	epT.I.P.S.® Reloads, BioBased, Biopur, 55 mm, orange	20 to 300 µL,	5 trays of 96 tips (480 tips)
041-033357	epT.I.P.S.® Reloads, BioBased, Biopur, 71 mm, blue	50 to 1000 µL	5 trays of 96 tips (480 tips)
041-033358	epT.I.P.S.® Reloads, BioBased, Biopur, 76 mm, green	50 to 1250 µL	5 trays of 96 tips (480 tips)
041-033359	epT.I.P.S.® Reloads, BioBased, Biopur, L, 103 mm, dark green	50 to 1250 µL	5 trays of 96 tips (480 tips)
	BioBased tips - Dualfilter - PCR Clean and Sterile Reloads		
041-033363	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, M, 40 mm, medium grey	0.1 to 10 µL	10 trays of 96 tips (960 tips)
041-033364	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, L, 46 mm, light grey	0.5 to 20 µL	10 trays of 96 tips (960 tips)
041-033365	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, 53 mm, yellow	2 to 20 µL	10 trays of 96 tips (960 tips)
041-033366	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, 53 mm, yellow	2 to 100 µL	10 trays of 96 tips (960 tips)
041-033367	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, 55 mm, yellow	2 to 200 µL	10 trays of 96 tips (960 tips)
041-033368	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, 55 mm, orange	20 to 300 µL	10 trays of 96 tips (960 tips)
041-033369	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, 76 mm, blue	50 to 1000 µL	10 trays of 96 tips (960 tips)
041-033370	ep Dualfilter T.I.P.S.® Reloads, BioBased, PCR clean and Sterile, L, 103 mm, dark green	50 to 1250 µL	10 trays of 96 tips (960 tips)
	BioBased tips - DualFilter - SealMax reloads		
041-033371	ep Dualfilter T.I.P.S.® SealMax, Reloads, BioBased, Biopur, L, 46 mm, light grey	0.5 to 20 µL	10 trays of 96 tips (960 tips)
041-033372	ep Dualfilter T.I.P.S.® SealMax, Reloads, BioBased, Biopur, 53 mm, yellow	2 to 100 µL	10 trays of 96 tips (960 tips)
041-033373	ep Dualfilter T.I.P.S.® SealMax, Reloads, BioBased, Biopur, 55 mm, yellow	2 to 200 µL	10 trays of 96 tips (960 tips)
041-033374	ep Dualfilter T.I.P.S.® SealMax, Reloads, BioBased, Biopur, 55 mm, orange	20 to 300 µL	10 trays of 96 tips (960 tips)
041-033375	ep Dualfilter T.I.P.S.® SealMax, Reloads, BioBased, Biopur, 76 mm, blue	50 to 1000 µL	10 trays of 96 tips (960 tips)
	BioBased tips - bulk		
041-033448	EpT.I.P.S.® tip 0.1-10 uL bulk BioBased dark grey x500x2	0.1 - 10 µL	Bulk : 1000 tips
041-033449	EpT.I.P.S.® tip 0.1-20 uL bulk BioBased light grey x500x2	0.1 - 20 µL	Bulk : 1000 tips
041-033450	EpT.I.P.S.® tip 0.5-20 uL, L, bulk BioBased light grey x500x2	0.5 - 20 µL L	Bulk : 1000 tips
041-033451	EpT.I.P.S.® tip 2-200 uL bulk BioBased yellow x500x2	20 - 300 µL	Bulk : 1000 tips
041-033452	EpT.I.P.S.® tip 20-300 uL bulk BioBased orange x500x2	2 - 200 µL	Bulk : 1000 tips
041-033453	EpT.I.P.S.® tip 50-1000 uL bulk BioBased blue x500x2	50 - 1000 µL	Bulk : 1000 tips
041-033454	EpT.I.P.S.® tip 50-1250 uL bulk BioBased green x500x2	50 - 1250 µL	Bulk : 1000 tips
041-033455	EpT.I.P.S.® tip 50-1250 uL, L, bulk BioBased dark green x250x4	50 - 1250 µL L	Bulk : 1000 tips
041-033456	EpT.I.P.S.® tip 250-2500 uL bulk BioBased red x100x5	0.25 - 2.5 mL	Bulk : 500 tips
041-033457	EpT.I.P.S.® tip 100-5000 uL bulk BioBased violet x100x5	0.1 - 5 mL	Bulk : 500 tips
041-033458	EpT.I.P.S.® tip 200-5000 uL, L, bulk BioBased violet x100x3	0.2-5 mL L	Bulk : 300 tips

THESE PIPETTE TIPS ARE MADE FROM AT LEAST 90% BIOBASED POLYPROPYLENE.

ACT. The Environmental Impact Factor Label

Excellent quality

- General quality certificates: compliance with general quality standards and requirements as well as product purity certificates
- ACT certificates: independent certification of the environmental impact factor awarded by My Green Lab

Support us in our global sustainability efforts:

- Collect the cardboard packaging from your tubes.
- Dispose of the packaging in the collection bins indicated by your organisation and region.







milian DUTSCHER GROUP

Chemin du Champ-des-Filles n°4A
Campus Spark
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
www.milian.com