milian®
DUTSCHER GROUP

BAMBANKER
Cryopreservation
Freezing medium

NIPPON Genetics EUROPE
INNOVATION FOR YOU

Cryopreservation
Bambanker freezing medium

Cell freezing medium, ready to use without SVF, with or without DMSO. Fast and long-term freezing at - 80 °C and product storage for 2 years at + 4 °C.

4 versions:

- Standard: for sensitive cells such as embryonic stem cells and primary cells
- Direct: no need for centrifugation, suitable for hybridomas and HTP
- HRM: contains human serum albumin, no animal-derived components, for research use only
- Without DMSO: for cell lines sensitive to DMSO, ready to use without dilution, direct freezing at -80 °C

Volume (ml)	100	20	20	10
Units/box	1	5	1	1
Standard version	523302B	523303	523304	523305
CHF	N/A -	N/A -	N/A -	N/A -
Direct version	-	-	523306	-
CHF	-	-	N/A -	-
HMR version	523631B	-	523308	-
CHF	N/A -	-	N/A -	-
DMSO-free version	-	-	523374	-
CHF	-	-	N/A -	-



523302B



523306



523308



523374

CONTACT YOUR LIFE SCIENCES EXPERT!



Deborah FEINTUCH
Life Sciences Specialist in CHR
dfeintuch@milian.com
+41 79 526 13 29



Matthias GROSS
Life Sciences Specialist in CHD
mgross@milian.com
+41 79 505 56 15