milian®
DUTSCHER GROUP

TIPS FOR LTS MICROPIPETTES


ClearLine®

ClearLine® LTS Compatible Pipet Tips

The new ClearLine® LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch™ pipettors, all with the quality you have come to expect from ClearLine®.

- **Comfort you can feel** : LTS Compatible Tips reduce the force required to seal and eject, enabling a streamlined experience from aspiration to dispensing. Avoid slamming with excessive force on the pipettor for a proper seal, eliminate the struggle to eject tips, and minimize hand and wrist strain after extensive pipetting
- **Convenience of choice** : LTS Compatible Tips are available in your favorite packaging formats, including racked, racked/filtered, tip station, reload stacks, and our new more sustainable solution: ClearLine® HybridRack™ Pipet Tips. If you're not using an LTS pipettor, no problem. We provide a complete range of universal pipet tips to choose from
- **Quality** : Our United States manufacturing facility is audited by third parties to ensure adherence to regulated manufacturing standards, including ISO 13485 and ISO 14001. Sterile products feature a Sterility Assurance Level (SAL) of 10^{-3}, and all products are made using USP Class VI resin
- **Peace of mind** : Our LTS Compatible pipet tips are RNase-/DNase-free, Human gDNA-free, PCR inhibition-free, and nonpyrogenic so you can pipet stress-free, knowing your samples are safe. Additionally, the tips are autoclavable at 121°C for 15 minutes

Effortless ejection, superior sealing, ergonomic excellence



milian®
DUTSCHER GROUP

Campus Spark
Chemin du Champ-des-Filles n°4A
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
milian@milian.com / www.milian.com

milian®
DUTSCHER GROUP

TIPS FOR LTS MICROPIPETTES

ClearLine®

ClearLine® filter and non-filter tips for LTS micropipettes

- Polypropylene filter and non-filter tips
- Compatible with Rainin LTS micropipettes
- 4 racks **to choose from**

DNase RNase · ADN humain · Inhibiteur de PCR · Pyro-gènes



ClearLine® Plus rack:

+ Recyclable plastic hinged rack* ;
+ colour-coded tray ;
+ rechargeable ;
+ stackable and autoclavable.




ClearLine® Plus refill :

+ Recyclable, colour-coded refill tray;
+ Refill for Clearline Plus and HybridRack racks.



TipStation:

+ 5 layers of 192 tips, minimum footprint ;
+ simply remove one layer after using the tips;
+ very stable basis ;
+ ideal for multichannel micropipettes ;
+ in recyclable polypropylene*.



HybridRack :

+ Rack and covering lid made from easily recyclable recycled cardboard ;
+ coloured tray in recyclable plastic* ;
+ 77 to 86% less plastic waste than a conventional plastic rack ;
+ Tray compatible with ClearlinePlus racks.

* see our recycling programme on milian.com

Non-filter Tips

						ClearLine® Plus rack				ClearLine® Plus refill		TipStation:				
						Sterile		Non-sterile		Non-sterile			Sterile		Non-sterile	
Photo number**	Vol. (µl)	Designation	Tray colour	Graduation (µl)	Units/ box	Part no.	CHF/ box	Part no.	CHF/ box	Part no.	CHF/ box	Tray colour	Part no.	CHF/ box	Part no.	CHF/ box
1	20	Fine tip, length 49.5 mm	Red	2, 10	10 x 96	**258380**	**NC -**	**258383**	**NC -**	**258386**	**NC -**	Blue	**258389**	**NC -**	**258391**	**NC -**
2	200	Tip length 50.2 mm	Yellow	10, 50, 100	10 x 96	**258381**	**NC -**	**258384**	**NC -**	**258387**	**NC -**	Blue	**258390**	**NC -**	**258392**	**NC -**
3	1000	Tip length 98.6 mm	Blue	200, 500, 1000	8 x 96	**258382**	**NC -**	**258385**	**NC -**	**258388**	**NC -**	-	-	-	-	-

** Refer to the picture of the model with filter (below) without taking this one into account.

Filter tips

Tips not shown to scale.	Num. photo	Vol. (µl)	Designation	Tray colour	Graduation (µl)	Units/ box	ClearLine® Plus rack Sterile		Hybridrack Sterile	
							Part no.	CHF/ box	Part no.	CHF/ box
	1	20	Fine tip, length 49.5 mm	Red	2, 10	10 x 96	**258377**	**NC -**	**258393**	**NC -**
	2	200	Tip length 50.2 mm	Yellow	10, 50, 100	10 x 96	**258378**	**NC -**	**258394**	**NC -**
	3	1000	Tip length 98.6 mm	Blue	200, 500, 1000	8 x 96	**258379**	**NC -**	**258395**	**NC -**

milian®
DUTSCHER GROUP

Campus Spark
Chemin du Champ-des-Filles n°4A
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
milian@milian.com / www.milian.com