

MIKRO 120



MIKRO 200 MIKRO 200 R



MIKRO 220 MIKRO 220 R



EBA 20



EBA 20S



EBA 21



EBA 270

| | | Mikro Zentrifugen | | | | | Kleine Zentrifugen | | | |
|---|---|---|---|---|---|---|---|---|---|---|
| | | MIKRO 120 | MIKRO 200 | MIKRO 200 R Gekühlt | MIKRO 220 | MIKRO 220 R Gekühlt | EBA 20 | EBA 20S | EBA 21 | EBA 270 |
| Winkelrotor | RZB max. (x g) | 17.530 | 21.382 | 21.382 | 31.514 | 31.514 | 3.461 | 6.153 | 25.718 | - |
| | U/min max. | 14.000 | 15.000 | 15.000 | 18.000 | 18.000 | 6.000 | 8.000 | 18.000 | - |
| Rotor frei | RZB max. (x g) | - | - | - | 18.516 | 18.516 | - | - | 2.999 | 2.254 |
| | U/min max. | - | - | - | 13.000 | 13.000 | - | - | 5.000 | 4.000 |
| Rotor mit Trommel | RZB max. (x g) | - | - | - | 14.171 | 14.171 | - | - | - | - |
| | U/min max. | - | - | - | 13.000 | 13.000 | - | - | - | - |
| Röhrchen Kolben Platten | Mikroröhrchen 0,2 - 0,5 ml (Ø 6 - 8 mm, Höhe 18 - 30 mm) | 24 | 30 | 30 | 192 | 192 | - | - | 30 | - |
| | Mikroröhrchen 1,5 ml (Ø 11 mm, Höhe 38 mm) | 24 | 30 | 30 | 60 | 60 | - | - | 30 | - |
| | Mikroröhrchen 2 ml (Ø 11 mm, Höhe 38 mm) | 24 | 30 | 30 | 60 | 60 | - | - | 30 | - |
| | Röhrchen mit rundem Boden 15 ml (Ø 17 mm, Höhe100 mm) | - | - | - | - | 12 | 8 | 8 | 12 | 6 |
| | Röhrchen mit rundem Boden 50 ml (Ø 34 mm, Höhe 100 mm) | - | - | - | - | 6 | - | - | 6 | - |
| | Röhrchen mit rundem Boden 100 ml | - | - | - | - | - | - | - | - | - |
| | Röhrchen mit rundem Boden 250 ml (Ø 65 mm, Höhe 115 mm) | - | - | - | - | - | - | - | - | - |
| | Blutröhrchen (Ø 13 mm, Höhe 75 / 100 mm) | - | - | - | - | 12 | 8 | 8 | 12 | 6 |
| | Blutröhrchen (Ø 16 mm, Höhe 75 / 100 mm) | - | - | - | - | 8/12 | -/ 8 | -/ 8 | 8 / 12 | 6 |
| | Konische Röhrchen 15 ml (Ø 17 mm, Höhe 120 mm) | - | - | - | - | 6 | 4 | 4 | 6 | - |
| | Konische Röhrchen 50 ml ohne Stehrand / mit Stehrand (Ø 29 mm, Höhe 115 mm) | - | - | - | - | 3 / - | - | - | 3 / - | - |
| | Konsiche Röhrchen breit 250 / 500 ml | - | - | - | - | - | - | - | - | - |
| | Grossvolumige Röhrchen 250 ml (Ø 61 mm, Höhe 122 mm) | - | - | - | - | - | - | - | - | - |
| | Grossvolumige Röhrchen 600 / 750 ml (Ø 93 mm, Höhe 124 mm / Ø 97 mm, Höhe 152 mm) | - | - | - | - | - | - | - | - | - |
| | Grossvolumige Röhrchen 1000 ml (Ø 96 mm, Höhe 176 mm) | - | - | - | - | - | - | - | - | - |
| | Grossvolumige Röhrchen 2000 ml (180 x 150 x 100 mm) | - | - | - | - | - | - | - | - | - |
| | Röhrchen für Blutbeutel 450 ml | - | - | - | - | - | - | - | - | - |
| | PCR-Strip 8 x 0,2 ml | - | 4 x 8 | 4 x 8 | 6 x 8 | 6 x 8 | - | - | - | - |
| | MTP-Platten (86 x 128 x 15 mm) | - | - | - | - | - | - | - | - | - |
| | Zytologie | - | - | - | - | - | - | - | - | - |



ROTOFIX 32 A



ROTINA 380 ROTINA 380 R



ROTINA 420 ROTINA 420 R



ROTANTA 460 RC



ROTIXA 50 S ROTIXA 50 RS



ROTO SILENTA 630 RS



UNIVERSAL 320 UNIVERSAL 320 R



ROTANTA 460



ROTANTA 460 R



ROTANTA 460 RF

| | | Labortisch-Zentrifugen | | | | | | | | | | | Zentrifugen auf Standfüssen | | |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | | ROTOFIX 32 A | UNIVERSAL 320 | UNIVERSAL 320 R | ROTINA 380 | ROTINA 380 R | ROTINA 420 | ROTINA 420 R | ROTANTA | ROTANTA 460 R | ROTANTA 460 RC | ROTANTA 460 RF | ROTIXA 50 S | ROTIXA 50 RS | ROTO SILENTA 630 RS |
| | | | | Gekühlt | | Gekühlt | | Gekühlt | 460 | Gekühlt | Gekühlt | Gekühlt | | Gekühlt | - |
| Winkelrotor | RZB max. (x g) | 4.226 | 21.382 | 21.382 | 24.400 | 24.400 | 24.400 | 24.400 | 24.400 | 24.400 | 24.400 | 24.400 | - | - | - |
| | U/min max. | 6.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | 15.000 | - | - | Gekühlt |
| Rotor frei | RZB max. (x g) | 2.701 | 4.193 | 4.193 | 4.863 | 4.863 | 4.166 | 4.166 | 5.858 | 6.446 | 6.446 | 6.446 | 5.530 | 5.530 | 6.498 |
| | U/min max. | 4.000 | 5.000 | 5.000 | 5.000 | 5.000 | 4.500 | 4.500 | 5.900 | 6.200 | 6.200 | 6.200 | 4.900 | 4.900 | 4.500 |
| Rotor mit Trommel | RZB max. (x g) | - | - | - | 18.866 | 18.866 | - | - | - | - | - | - | - | - | - |
| | U/min max. | - | - | - | 15.000 | 15.000 | - | - | - | - | - | - | - | - | - |
| Röhrchen Kolben Platten | Mikroröhrchen 0,2 - 0,5 ml (Ø 6 - 8 mm, Höhe 18 - 30 mm) | - | 30 | 30 | 192 | 192 | 30 | 30 | 30 | 30 | 30 | 30 | - | - | - |
| | Mikroröhrchen 1,5 ml (Ø 11 mm, Höhe 38 mm) | 36 | 48 | 48 | 160 | 160 | 160 | 160 | 240 | 240 | 240 | 240 | 336 | 336 | - |
| | Mikroröhrchen 2 ml (Ø 11 mm, Höhe 38 mm) | 36 | 48 | 48 | 160 | 160 | 160 | 160 | 240 | 240 | 240 | 240 | 336 | 336 | - |
| | Röhrchen mit rundem Boden 15 ml (Ø 17 mm, Höhe 100 mm) | 32 | 32 | 32 | 52 | 52 | 72 | 72 | 148 | 148 | 148 | 148 | 120 | 120 | 180 |
| | Röhrchen mit rundem Boden 50 ml (Ø 34 mm, Höhe 100 mm) | 6 | 6 | 6 | 12 | 12 | 16 | 16 | 16 | 16 | 16 | 16 | 24 | 24 | 36 |
| | Röhrchen mit rundem Boden 100 ml | 4 | 4 | 4 | 4 | 4 | 12 | 12 | 12 | 12 | 12 | 12 | 16 | 16 | 24 |
| | Röhrchen mit rundem Boden 250 ml (Ø 65 mm, Höhe 115 mm) | - | - | - | 4 | 4 | 4 | 4 | 6 | 6 | 6 | 6 | 6 | 6 | 6 |
| | Blutröhrchen (Ø 13 mm, Höhe 75 / 100 mm) | 32 / 20 | 32 | 32 | 64 | 64 | 84 / 104 | 84 / 104 | 196 | 196 | 196 | 196 | 120 | 120 | 180 |
| | Blutröhrchen (Ø 16 mm, Höhe 75 / 100 mm) | 28 / 16 | 28 | 28 | 52 | 52 | 72 | 72 | 148 | 148 | 148 | 148 | 120 | 120 | 180 |
| | Konische Röhrchen 15 ml (Ø 17 mm, Höhe 120 mm) | 32 | 32 | 32 | 36 | 36 | 52 | 52 | 96 | 96 | 96 | 96 | 92 | 92 | 138 |
| | Konische Röhrchen 50 ml ohne Stehrand / mit Stehrand (Ø xH : 29 x 115 mm) | 8 | 8 | 8 | 16 | 16 | 20 | 20 | 40 | 40 | 40 | 40 | 32 / 24 | 32 / 24 | 48 / 36 |
| | Konische Röhrchen breit 250 / 500 ml | - | - | - | - | - | - | - | 4 | 4 | 4 | 4 | 4 | 4 | 6 |
| | Grossvolumige Röhrchen 250 ml (Ø 61 mm, Ht 122 mm) | - | - | - | 4 | 4 | 4 | 4 | 8 | 8 | 8 | 8 | 6 | 6 | 6 |
| | Grossvolumige Röhrchen 600 / 750 ml (Ø x H : 93 x 124 mm / Ø x H : 97 x 152 mm) | - | - | - | - | - | 4 x 600 | 4 x 600 | 4 x 750 | 4 x 750 | 4 x 750 | 4 x 750 | 4 | 4 | 6 |
| | Grossvolumige Röhrchen 1000 ml (Ø x H: 96 x 176 mm) | - | - | - | - | - | - | - | - | - | - | - | 4 | 4 | 6 |
| | Grossvolumige Röhrchen 2000 ml (180 x 150 x 100 mm) | - | - | - | - | - | - | - | - | - | - | - | - | - | 6 |
| | öhrchen für Blutbeutel 450 ml | - | - | - | - | - | - | - | 4 | 4 | 4 | 4 | 4 | 4 | 12 |
| | PCR-Strip 8 x 0,2 ml | - | 24 x 8 | 24 x 8 | 24 x 8 | 24 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | 48 x 8 | - |
| | MTP-Platten (86 x 128 x 15 mm) | - | 10 | 10 | 10 | 10 | 16 | 16 | 24 | 24 | 24 | 24 | 16 | 16 | - |
| | Zytologie | 8 | 8 | 8 | 8 | 8 | 8 | 8 | 16 | 16 | 16 | 16 | 12 | 12 | - |